

# Accela Open Autosampler

## User Guide for LC Devices

60357-97101 Revision D December 2013



**Thermo**  
SCIENTIFIC

© 2013 Thermo Fisher Scientific Inc. All rights reserved.

Exactive, Foundation, LCQuan, and MSQ are trademarks, and Accela, LCQ, LCQ Fleet, LTQ, LXQ, Thermo Scientific, and Xcalibur are registered trademarks of Thermo Fisher Scientific Inc. in the United States.

The following are registered trademarks in the United States and other countries: Adobe and Reader are registered trademarks of Adobe Systems Incorporated. Agilent, ChemStation, and Instant Pilot are registered trademarks of Agilent Technologies Inc. Empower and MassLynx are registered trademarks of Waters Technologies Corporation. Microsoft and Windows are registered trademarks of Microsoft Corporation.

The following are registered trademarks in the United States and possibly other countries: Analyst is a registered trademark of AB Sciex Pte. Ltd. Cheminert is a registered trademark of Valco Instruments Company, Inc. Kalrez is a registered trademark of E.I. du Pont de Nemours & Company. Ryton is a registered trademark of Chevron Phillips Chemical Company. Simriz is a registered trademark of Carl Freudenberg KG Corporation.

PAL is a trademark of CTC Analytics AG. PEEK is a trademark of Victrex plc.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo Fisher Scientific Inc. provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Thermo Fisher Scientific Inc.

The contents of this document are subject to change without notice. All technical information in this document is for reference purposes only. System configurations and specifications in this document supersede all previous information received by the purchaser.

This document is not part of any sales contract between Thermo Fisher Scientific Inc. and a purchaser. This document shall in no way govern or modify any Terms and Conditions of Sale, which Terms and Conditions of Sale shall govern all conflicting information between the two documents.

Release history: Revision A, Apr 2010; Revision B, Mar 2011; Revision C, Jan 2013; Revision D, Dec 2013

Software version: Thermo Foundation (see page xi); Thermo Xcalibur 2.1 or later; Thermo LC Devices 2.5.0 or later. Firmware version (see page xii).

**For Research Use Only. Not for use in diagnostic procedures.**

## Regulatory Compliance

Thermo Fisher Scientific performs complete testing and evaluation of its products to ensure full compliance with applicable domestic and international regulations. When the system is delivered to you, it meets all pertinent electromagnetic compatibility (EMC) and safety standards as described in the next section or sections by product name.

Changes that you make to your system may void compliance with one or more of these EMC and safety standards. Changes to your system include replacing a part or adding components, options, or peripherals not specifically authorized and qualified by Thermo Fisher Scientific. To ensure continued compliance with EMC and safety standards, replacement parts and additional components, options, and peripherals must be ordered from Thermo Fisher Scientific or one of its authorized representatives.

### EMC Directive 2004/108/EC

EMC compliance has been evaluated by CTC Analytics for the HTC PAL™ autosampler.

IEC 61326-1:2005

IEC 61326-2-6:2005

EN 61326-1:1997, A1:1998

CISPR 22:2005, A1:2005, A2:2006

FCC Class A, CFR 47 Part 15:2003

### Low Voltage Safety Compliance

Low Voltage Safety Compliance has been evaluated by CTC Analytics for the HTC PAL autosampler.

This device complies with Low Voltage Directive 2006/95/EC and harmonized standard EN 61010-1:2001, IEC 61010-1:2001, ANSI/UL 61010 A-1:2004, CAN/CSA 22.2 61010-1:2004.

## FCC Compliance Statement

THIS DEVICE COMPLIES WITH PART 15 OF THE FCC RULES. OPERATION IS SUBJECT TO THE FOLLOWING TWO CONDITIONS: (1) THIS DEVICE MAY NOT CAUSE HARMFUL INTERFERENCE, AND (2) THIS DEVICE MUST ACCEPT ANY INTERFERENCE RECEIVED, INCLUDING INTERFERENCE THAT MAY CAUSE UNDESIRE OPERATION.



**CAUTION** Read and understand the various precautionary notes, signs, and symbols contained inside this manual pertaining to the safe use and operation of this product before using the device.



## Notice on Lifting and Handling of Thermo Scientific Instruments

For your safety, and in compliance with international regulations, the physical handling of this Thermo Fisher Scientific instrument *requires a team effort* to lift and/or move the instrument. This instrument is too heavy and/or bulky for one person alone to handle safely.

## Notice on the Proper Use of Thermo Scientific Instruments

In compliance with international regulations: This instrument must be used in the manner specified by Thermo Fisher Scientific to ensure protections provided by the instrument are not impaired. Deviations from specified instructions on the proper use of the instrument include changes to the system and part replacement. Accordingly, order replacement parts from Thermo Fisher Scientific or one of its authorized representatives.

## Notice on the Susceptibility to Electromagnetic Transmissions

Your instrument is designed to work in a controlled electromagnetic environment. Do not use radio frequency transmitters, such as mobile phones, in close proximity to the instrument.

For manufacturing location, see the label on the instrument.

## WEEE Compliance

This product complies with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2002/96/EC. It is marked with the following symbol:



Thermo Fisher Scientific has contracted with one or more recycling or disposal companies in each European Union (EU) Member State, and these companies should dispose of or recycle this product. For further information, go to [www.thermoscientific.fr/rohs](http://www.thermoscientific.fr/rohs).

## WEEE Konformität

Dieses Produkt entspricht der EU Waste Electrical & Electronic Equipment (WEEE) Richtlinie 2002/96/EC. Es ist mit dem folgenden Symbol gekennzeichnet:



Thermo Fisher Scientific hat Vereinbarungen mit Verwertungs-/Entsorgungsfirmen in allen EU-Mitgliedsstaaten getroffen, damit dieses Produkt durch diese Firmen wiederverwertet oder entsorgt werden kann. Weitere Informationen finden Sie unter [www.thermoscientific.de/rohs](http://www.thermoscientific.de/rohs).

## Conformité DEEE

Ce produit est conforme avec la directive européenne (2002/96/EC) des Déchets d'Equipements Electriques et Electroniques (DEEE). Il est marqué par le symbole suivant:



Thermo Fisher Scientific s'est associé avec une ou plusieurs sociétés de recyclage dans chaque état membre de l'Union Européenne et ce produit devrait être collecté ou recyclé par celle(s)-ci. Pour davantage d'informations, rendez-vous sur la page [www.thermoscientific.fr/rohs](http://www.thermoscientific.fr/rohs).

# Contents

	<b>Preface</b> .....	<b>ix</b>
	Related Documentation .....	x
	System Requirements .....	xi
	Supported Firmware Versions .....	xii
	Cautions and Special Notices .....	xiii
	Contacting Us .....	xiii
<b>Chapter 1</b>	<b>Getting Connected</b> .....	<b>1</b>
	Ordering Information .....	1
	Setting Up an LC/MS System with the Autosampler .....	2
	Connecting the Accela Open Autosampler to the Data System Computer .....	4
	System Synchronization Connections .....	4
	Connecting the Accela Open Autosampler Interconnect Cable .....	6
	Connecting the Adapter Cable for an Accela Detector or TSQ Series Mass Spectrometer .....	8
	Connecting the Accela Pump Adapter Cable .....	12
<b>Chapter 2</b>	<b>Getting Started</b> .....	<b>15</b>
	Configuring the Accela Open Autosampler .....	15
	Using the Autosampler Control Terminal .....	19
	Selecting the Tray Type .....	20
	Changing the Pulse Time Setting .....	21
	Specifying the Instrument Method Parameters .....	22
	Viewing the Method Summary .....	29
<b>Chapter 3</b>	<b>Creating Custom Templates and Macros</b> .....	<b>31</b>
	Using the Template Editor to Create Custom Templates .....	31
	Testing a Custom Template .....	34
	Viewing the Template Summary .....	37
	Standard Macros .....	39
	Using the Macro Editor to Create Custom Macros .....	39
	Defining Variables .....	44

<b>Chapter 4</b>	<b>Using the Direct Controls to Operate the Autosampler</b> . . . . .	<b>47</b>
	Opening the Direct Control Dialog Box . . . . .	47
	Making Single Injections from the Tune Window . . . . .	48
	Changing the Tray Configuration . . . . .	50
	Specifying the Location of the Direct Control Methods . . . . .	52
	Wrapping Direct Control Injections with Additional Macros . . . . .	53
<b>Chapter 5</b>	<b>Using Dynamic Load and Wash (DLW)</b> . . . . .	<b>55</b>
	Overview . . . . .	55
	Cycle Composer Macros or ICC Cycles . . . . .	57
	General Considerations . . . . .	58
	Priming the Solvent Lines, Wash1 and Wash2 . . . . .	59
	Standard DLW Injection Cycle . . . . .	60
	Fast DLW Injection Cycle . . . . .	63
	Operating Dynamic Load and Wash (DLW) . . . . .	64
	Priming the Solvent Lines . . . . .	65
	Location of Solvent and Waste Bottles . . . . .	65
	DLW Pumps . . . . .	66
	DLW Actuator/Solenoid . . . . .	66
	DLW Cycle Step-by-Step . . . . .	68
	Cycle for Standard Injection . . . . .	68
	Additional Valve Toggle Step to DLW Standard Cycle . . . . .	82
	Cycle for Fast Injection . . . . .	90
<b>Appendix A</b>	<b>Sample Trays</b> . . . . .	<b>103</b>
	Tray Types . . . . .	104
	Tray Position . . . . .	106
	Sample Position . . . . .	107
	<b>Index</b> . . . . .	<b>109</b>

# Preface

This guide describes how to connect a Thermo Scientific™ Accela™ Open Autosampler to other devices in the liquid chromatography/mass spectrometry (LC/MS) system and how to control the autosampler from a Thermo Scientific MS data system, such as Thermo Xcalibur™.

The Accela Open Autosampler manufactured by CTC Analytics is supplied by Thermo Fisher Scientific as part of a solution for high-throughput LC/MS applications.

## Contents

- [Related Documentation](#)
- [System Requirements](#)
- [Supported Firmware Versions](#)
- [Cautions and Special Notices](#)
- [Contacting Us](#)

### ❖ To suggest changes to the documentation or to the Help

Complete a brief survey about this document by clicking the button below.  
Thank you in advance for your help.



## Related Documentation

Thermo LC Devices includes Help and these manuals as PDF files:

- *Accela Open Autosampler Hardware Manual*
- *Accela Open Autosampler User Guide*

### ❖ To view product manuals

From the Microsoft™ Windows™ taskbar, do the following:

- For an LC instrument controlled by a Thermo software application, choose **Start > All Programs > Thermo Instruments > Manuals > LC Devices > Accela.**
- For the Xcalibur manual set, choose **Start > All Programs > Thermo Xcalibur > Manuals > Xcalibur.**

For access to the application Help, follow this procedure.

### ❖ To view application-specific Help

- From the Thermo Xcalibur Instrument Setup window, choose **Help > Accela Open AS Help.**
- If information about setting parameters is available for a specific view, page, or dialog box, click **Help** or press the F1 key for information about setting parameters.

For more information, visit [www.thermoscientific.com](http://www.thermoscientific.com).

## System Requirements

Ensure that the system meets these minimum requirements.

**IMPORTANT** Before you install the device driver, ensure that the data system computer has a compatible version of the Thermo Foundation™ platform as noted in the *Thermo LC Devices x.x.x Release Notes*.

System	Minimum requirements										
Computer	<ul style="list-style-type: none"> <li>• 2 GHz processor with 1 GB RAM</li> <li>• DVD drive</li> <li>• 80 GB available on drive C</li> <li>• Video card and monitor capable of 1280 × 1024 resolution</li> <li>• NTFS format</li> </ul>										
Software	<ul style="list-style-type: none"> <li>• Adobe™ Reader™ 9.0</li> <li>• Microsoft Windows operating system:               <ul style="list-style-type: none"> <li>– Windows 7 Professional (32-bit and 64-bit<sup>a</sup>)</li> <li>– Windows XP Workstation SP3</li> </ul> </li> <li>• Thermo Scientific software:               <table border="1" data-bbox="678 1058 1240 1293"> <thead> <tr> <th>LC Devices</th> <th>Foundation</th> </tr> </thead> <tbody> <tr> <td>2.6.0 and later</td> <td>1.0.2 SP2 and later</td> </tr> <tr> <td>2.5.0 SP3</td> <td>1.0.2 SP2</td> </tr> <tr> <td>2.5.0 SP1 or SP2</td> <td>2.0</td> </tr> <tr> <td>2.5.0</td> <td>1.0.2 SP2</td> </tr> </tbody> </table> </li> </ul>	LC Devices	Foundation	2.6.0 and later	1.0.2 SP2 and later	2.5.0 SP3	1.0.2 SP2	2.5.0 SP1 or SP2	2.0	2.5.0	1.0.2 SP2
LC Devices	Foundation										
2.6.0 and later	1.0.2 SP2 and later										
2.5.0 SP3	1.0.2 SP2										
2.5.0 SP1 or SP2	2.0										
2.5.0	1.0.2 SP2										

<sup>a</sup>LC Devices 2.8.0 and later are compatible with Windows 7, 32-bit and 64-bit.


(Windows 7 only) If you receive a server failure error when you try to open the Xcalibur Instrument Setup window, follow the next procedure.

### ❖ To resolve a server failure for the Xcalibur data system

1. Verify that the installed versions of Foundation platform and LC Devices are compatible (see [System Requirements](#)).
2. If the installed LC Devices software is compatible with Foundation, go to [step 3](#). If it is not compatible, do the following:
  - a. Use the Windows Control Panel to uninstall all of the modules from LC Devices.
  - b. Install the compatible version of LC Devices.
  - c. Restart the data system computer.

3. If the installed LC Devices is compatible with Foundation, do the following:
  - a. Open the Windows Control Panel.
  - b. In the top Search box, type **Change User Account Control Settings**, and then select this link to open the User Account Control Settings dialog box.
  - c. Move the slider to the **Default** position (Figure 1).
  - d. Click **OK**.

**Figure 1.** User Account Control Settings dialog box



## Supported Firmware Versions

LC Devices supports firmware version 4.1.2 for the Thermo Scientific Accela Open Autosampler.

## Cautions and Special Notices

Make sure you follow the cautions and special notices presented in this guide. Cautions and special notices appear in boxes; those concerning safety or possible system damage also have corresponding caution symbols.

This guide uses the following types of cautions and special notices.



**CAUTION** Highlights hazards to humans, property, or the environment. Each CAUTION notice is accompanied by an appropriate CAUTION symbol.

**IMPORTANT** Highlights information necessary to prevent damage to software, loss of data, or invalid test results; or might contain information that is critical for optimal performance of the system.

**Note** Highlights information of general interest.

**Tip** Highlights helpful information that can make a task easier.

## Contacting Us

There are several ways to contact Thermo Fisher Scientific for the information you need.

### ❖ To contact Technical Support

Phone	800-532-4752
Fax	561-688-8736
E-mail	<a href="mailto:us.techsupport.analyze@thermofisher.com">us.techsupport.analyze@thermofisher.com</a>

### ❖ To contact Customer Service for ordering information

Phone	800-532-4752
Fax	561-688-8731
E-mail	<a href="mailto:us.customer-support.analyze@thermofisher.com">us.customer-support.analyze@thermofisher.com</a>
Web site	<a href="http://www.thermoscientific.com/en/support-landing/support.html">http://www.thermoscientific.com/en/support-landing/support.html</a>

❖ **To copy manuals from the Internet**

Go to [mssupport.thermo.com](https://mssupport.thermo.com), agree to the Terms and Conditions, and then click **Customer Manuals** in the left margin of the window.

❖ **To suggest changes to the Help**

- Fill out a reader survey online at [www.surveymonkey.com/s/PQM6P62](https://www.surveymonkey.com/s/PQM6P62).
- Send an e-mail message to the Technical Publications Editor at [techpubs-lcms@thermofisher.com](mailto:techpubs-lcms@thermofisher.com).

# Getting Connected

The Thermo Scientific Accela Open Autosampler is a high-throughput autosampler manufactured by CTC Analytics and supplied by Thermo Fisher Scientific. For high-throughput LC/MS applications, use the Accela Open Autosampler, a Thermo Scientific or an Agilent™ liquid chromatography (LC) pump, and a Thermo Scientific mass spectrometer (MS).

This chapter describes the contact closure connections between the Accela Open Autosampler, one or two Thermo Scientific LC pumps or an Agilent pump, a Thermo Scientific mass spectrometer, and the Xcalibur data system computer.



**CAUTION** Follow all of the recommendations given in the Safety Information section of the *CTC Analytics PAL System User Manual*. The Safety Information section includes information on the autosampler's electrical hazards, lithium battery, and safety labels.

## Contents

- [Ordering Information](#)
- [Setting Up an LC/MS System with the Autosampler](#)
- [Connecting the Accela Open Autosampler to the Data System Computer](#)
- [System Synchronization Connections](#)

## Ordering Information

Thermo Fisher Scientific supplies these cables with the Accela Open Autosampler:

- Serial communication cable that connects the autosampler to the data system computer
- Contact closure cable that connects the autosampler to the Thermo Scientific LC pump and Thermo Scientific mass spectrometer

To order the HTC PAL or HTS PAL autosampler from Thermo Fisher Scientific, see [Table 1](#).

**Table 1.** Accela Open Autosampler ordering information

Description	Part number
HTC PAL	OPTON 13009
HTC PAL with the 98 × 2 mL tray holder	OPTON 13010
HTS PAL	OPTON 13008

[Table 2](#) lists the contact closure cables used to connect an Agilent pump to a Thermo Scientific mass spectrometer and the autosampler. You can order the DB15 connector × 15 wire cable from Agilent or Thermo Fisher Scientific.

**Table 2.** Contact closure cables used to connect an Agilent pump

Description	Supplier	Part number
DB15 connector × 15 wire cable	Thermo Fisher Scientific	00012-27716
DB15 connector × 15 wire cable	Agilent	G1103-61611
PAL Interface Cable for APG Remote DB15 connector × DB9 connector cable	CTC Analytics	CBL 7890

## Setting Up an LC/MS System with the Autosampler

Typical stackable setups include placing the autosampler on top of the two pumps, in a dual-pump setup, and to the left of the mass spectrometer. Ensure that the stackable area for the autosampler is level, and that system cables behind the autosampler have adequate space.




**CAUTION** To prevent damage to the injection unit during operation, place the autosampler on a level surface.

For information on how to assemble the Accela Open Autosampler and connect the solvent lines, refer to the *Accela Open Autosampler Hardware Manual*.

[Figure 2](#) shows the front view of the Accela Open Autosampler. During installation, the Thermo Fisher Scientific field service engineer connects the following:

- LC system to the mass spectrometer
- Solvent lines from the injection valve to the LC pump and the Thermo Scientific mass spectrometer

**Figure 2.** Accela Open Autosampler




# Connecting the Accela Open Autosampler to the Data System Computer

The Accela Open Autosampler communicates with the data system computer through the supplied RS-232 serial communication cable.

### ❖ To connect the autosampler to the data system computer

Using the RS-232 serial cable supplied with the autosampler, connect the SER 1 port on the back of the autosampler to the serial port on the back of the data system computer (Figure 3).

**Figure 3.** Data system connection (firmware 4.1.2)



## System Synchronization Connections

The system interconnect cables and adapter cables that synchronize the run signals for an LC or LC/MS system with an Accela Open Autosampler depend on the Accela pump model, the mass spectrometer model, and whether the LC system includes an Accela detector.

Table 3 lists the system interconnect and adapter cables that Thermo Fisher Scientific supplies with the Accela Open Autosampler or by special order.

**Table 3.** Accela Open Autosampler contact closure and adapter cables

Cable	Function	Part number
Accela Open Autosampler interconnect cable <sup>a</sup>	Required for all hardware configurations	60157-63024
Accela detector adapter cable	Required for hardware configurations that include one or more of the following: <ul style="list-style-type: none"> <li>• Accela detector</li> <li>• TSQ Series or Exactive™ Series MS</li> </ul>	60157-63026
Accela Pump adapter cable <sup>b</sup>	Required for the discontinued Accela Pump	60157-63022

<sup>a</sup> Supplied with the autosampler

<sup>b</sup> Available by special order

**Note** The autosampler ships with an unlabeled cable that is equivalent to the Accela Open Autosampler interconnect cable (P/N 60157-63024). This interconnect cable has three connectors:

- The DB15 connector plugs into the INTERFACE 1 port on the back of the autosampler.
- The DB25 connector plugs into the D-Sub 25 port on the back of an Accela 600 or 1250 Pump.
- The 2-pin connector connects to the Accela detector adapter cable.


To connect the required cables listed in [Table 3](#), follow these procedures as appropriate:

- [Connecting the Accela Open Autosampler Interconnect Cable](#)
- [Connecting the Adapter Cable for an Accela Detector or TSQ Series Mass Spectrometer](#)
- [Connecting the Accela Pump Adapter Cable](#)

## Connecting the Accela Open Autosampler Interconnect Cable

Figure 4 shows the Accela Open Autosampler interconnect cable.


**Figure 4.** Accela Open Autosampler interconnect cable (P/N 60157-63024)



### ❖ To connect the Accela Open Autosampler interconnect cable to the LC or LC/MS system modules

1. Connect the DB15 connector (near label TO A/S) to the INTERFACE 1 port on the back panel of the autosampler (Figure 5).
2. Depending on the Accela pump model, do one of the following:
  - If your LC system includes an Accela 600 or 1250 Pump, connect the DB25 connector (near label TO PUMP) to the D-Sub25 port on the back of the pump (Figure 5).
  - If your LC system includes an Accela Pump, see “Connecting the Accela Pump Adapter Cable” on page 12.
3. Depending on the LC detector or MS setup, do one of the following:
  - If your LC or LC/MS system includes an Accela detector or a TSQ Series or Exactive MS, see “Connecting the Adapter Cable for an Accela Detector or TSQ Series Mass Spectrometer” on page 8.
  - If your LC/MS system includes an MSQ Plus™ MS or an LTQ™ Series MS, but does not include an Accela detector, go to step 4.
4. Connect the cable labeled TO MSQ/LTQ/LXQ to the MS. Do one of the following:
  - For the MSQ Plus MS, connect the cable to the User I/O Start In pins on the back of the MS (Figure 5).
  - For the LTQ Series MS, connect the cable to the Peripheral Control Start In pins on the right side of the MS (Figure 5).


**Figure 5.** Accela Open Autosampler interconnect cable connections (back view of the LC/MS modules)



## Connecting the Adapter Cable for an Accela Detector or TSQ Series Mass Spectrometer

Use the Accela detector and mass spectrometer adapter cable (Figure 6) to connect the Accela Open Autosampler to an Accela detector, a TSQ or Exactive mass spectrometer, or both.


**Figure 6.** Adapter cable for an Accela detector and a TSQ Series MS (P/N 60157-63026)



### ❖ To connect the LC detector and MS adapter cable

1. If you have not already done so, connect the Accela Open Autosampler interconnect cable to the autosampler and the Accela 600 or 1250 Pump (see [“Connecting the Accela Open Autosampler Interconnect Cable”](#) on page 6).
2. Connect the adapter cable to the interconnect cable by plugging the 2-pin male connector labeled **TO 60157-63024** on the adapter cable into the 2-pin female connector labeled **TO MSQ/LTQ/LXQ** on the interconnect cable (Figure 7).

**Figure 7.** Accela Open Autosampler interconnect cable to detector adapter cable connection




3. Using the adapter cable connectors, make the appropriate connections to an Accela detector, a Thermo Scientific MS, or both as follows:


- For an Accela detector (UV/Vis or PDA), connect the 8-pin connector attached to the cable labeled **TO PDA** to pins 1–8 on the back of the detector.
- For an MSQ Plus MS, connect the 2-pin connector attached to the cable labeled **TO LTQ/LXQ/MSQ** to the Start In pins on the back of the MS.
- For an LTQ Series MS, connect the 2-pin connector attached to the cables labeled **TO LTQ/LXQ/MSQ** and **TO LCQ FLEET** to the Start In pins on the right side of the MS.
- For a TSQ Series or Exactive Series MS, connect the cable labeled **TSQ/EXACTIVE** to the Start In pins on the right side of the MS.

Figure 8 shows the contact closure connections for a stand-alone Accela LC system, and Figure 9 on page 11 shows the contact closure connections for an LC/MS system.

**Figure 8.** Contact closure connections for a stand-alone Accela LC system with an Accela Open Autosampler



**Figure 9.** Accela Open Autosampler interconnect cable and detector adapter cable connections




## Connecting the Accela Pump Adapter Cable

Use the Accela Open Autosampler interconnect cable (Figure 4 on page 6) and the Accela Pump adapter cable (Figure 10) to interconnect the system modules for the following hardware configurations:

- Accela Open Autosampler, Accela Pump, and MSQ Plus mass spectrometer
- Accela Open Autosampler, Accela Pump, and a Thermo Scientific mass spectrometer


**Figure 10.** Accela Pump adapter cable



### ❖ To connect the Accela Pump adapter cable

1. If you have not already done so, connect the DB15 connector of the Accela Open Autosampler interconnect cable to the INTERFACE 1 port on the autosampler.
2. Connect the DB25 connector of the Accela Pump adapter cable to the DB25 connector of the Accela Open Autosampler interconnect cable (Figure 11).

**Figure 11.** Accela Open Autosampler Interconnect cable and Accela Pump adapter cable connection



3. Connect the 8-pin connector of the Accela Pump adapter cable to pins 1–8 on the back panel of the Accela Pump (Figure 12).
4. Connect the 2-pin connector of the Accela Open Autosampler interconnect cable to the MS (see “Connecting the Accela Open Autosampler Interconnect Cable” on page 6).





## Getting Started

This chapter describes how to add the autosampler to the Xcalibur instrument configuration, how to specify the available tray types, and how to specify the injection parameters for the autosampler from the Xcalibur data system.

### Contents

- [Configuring the Accela Open Autosampler](#)
- [Using the Autosampler Control Terminal](#)
- [Specifying the Instrument Method Parameters](#)

## Configuring the Accela Open Autosampler

The Accela Open Autosampler device driver provides control of the autosampler from the Xcalibur data system. To control the autosampler from the data system, add the autosampler to the Thermo Foundation Instrument Configuration window and specify its configuration options.


### ❖ To add the autosampler to the Foundation instrument configuration

1. Turn on the autosampler.
2. From the Microsoft Windows taskbar, choose **Start > All Programs > Thermo Foundation x.x > Instrument Configuration** to open the Thermo Foundation Instrument Configuration window ([Figure 13](#)).

## 2 Getting Started

### Configuring the Accela Open Autosampler

**Figure 13.** Thermo Foundation Instrument Configuration window




3. Under Available Devices, double-click the autosampler icon to add it to the Configured Devices area.

Use the Pal Configuration dialog box to select the PC Com port where the autosampler is connected, specify macros to import, and select sample trays or stacks.

#### ❖ To specify the configuration options for the autosampler

1. Follow the procedure [“To add the autosampler to the Foundation instrument configuration”](#) on [page 15](#).
2. In the Instrument Configuration window, under Configured Devices, double-click the autosampler icon, and then do the following on the displayed page ([Figure 14](#)).

**Figure 14.** Pal Configuration dialog box



3. Under Communications, do the following:

- a. In the Com Port list, select the communication port on the data system computer that the autosampler connects to. Or, click **Auto Select**.

**Note** The auto-select feature works only if the autosampler is connected to the data system computer and is turned on.

- b. If the autosampler is already connected to the data system computer, click **Check** to verify the connection.
4. Under Look Ahead Injections, to shorten the time between injections, select the **Enable** check box.

**Tip** Do not modify or delete a sequence while the data system is running with the Look Ahead Injections enabled. Use the stop function to finish the current injection before you modify or delete a sequence.

5. Under Method Template\Macro Editor, to make the Method Template and Macro Editor dialog boxes available in the Instrument Setup window, select the **Enable** check box.
6. Under Cycle Composer Macro Importer, import the macros that you plan to use to create injection methods as follows:
  - a. In the Import Path box, type the complete path to the location of the macros. Or, click **Browse** to locate the directory that contains the appropriate macros.
  - b. In the Overwrite Existing Macros list, select one of the following:
    - **Never**—The autosampler does not import the specified macro if it has the same file name as an existing macro.

- **Ask**—The autosampler asks you if you want to overwrite existing macros.
  - **Always**—Imported macros overwrite existing macros.
- c. Click **Import**.
7. Under Active Trays, select the trays that you want to use.

You can use trays that have a Yes in the Allow Use column. Clicking the tray name clears the Allow Use cell and removes the tray from the list of available trays in the Xcalibur data system.

**IMPORTANT** Use the autosampler control terminal to specify the tray information. You can change the tray information when you control the autosampler from the Xcalibur data system, but these changes are temporary. When you close and reopen the Xcalibur data system, the tray information resets to the autosampler control terminal settings.

8. Click **OK**.
9. Configure the other instruments in the LC/MS system, such as the LC pump and the mass spectrometer.
10. In the Thermo Foundation Instrument Configuration window, click **Done**.

Table 4 lists the parameters for the Pal Configuration dialog box.

**Table 4.** Pal Configuration dialog box parameters (Sheet 1 of 2)

Parameter	Description
<b>Communications</b>	
Com Port	Select the computer port where you plug the autosampler communication cable.
Auto Select	Instructs the system to automatically determine the computer port that connects to the autosampler through the autosampler communication cable.
Check	Instructs the system to verify the computer port that connects to the autosampler through the autosampler communication cable.
<b>Look Ahead Injections</b>	
Enable	Select this check box to enable Look Ahead Injections.  When enabled, the autosampler sets up for the next injection during the current run. The setup for the next injection involves carrying out the steps up to, but not including, the injection.
<b>Method Template \ Macro Editor</b>	

**Table 4.** Pal Configuration dialog box parameters (Sheet 2 of 2)

Parameter	Description
Enable	Select this check box when you want the Method Template and Macro Editor to be accessed through Instrument Setup. Ensuring the check box is clear prevents access to the Method Template and the Macro Editor.
<b>Cycle Composer Macro Importer</b>	
Import Path	Enter the path to the Cycle Macro to be imported.
Overwrite Existing Macros	Select Never, Ask, or Always.  The selection determines if existing macros are to be overwritten when a new macro is imported.
Import	Import the chosen macro.
Browse	Search through files and folders to find a macro.
<b>Active Trays</b>	
Tray Name	Hardware configured tray names.
Stack	Stack location of a given tray name.
Allow Use	Click to switch between Yes and blank. <ul style="list-style-type: none"> <li>• “Yes” means that the PAL system can access a given tray for obtaining samples (use is allowed).</li> <li>• A blank means that the PAL system cannot access the tray for obtaining samples (that is, use is not allowed).</li> </ul>

## Using the Autosampler Control Terminal

For general information about using the autosampler control terminal, refer to the *Accela Open Autosampler Hardware Manual*.

To select the sample trays and to change the default pulse time setting to 4.0 seconds, follow these procedures:


- [Selecting the Tray Type](#)
- [Changing the Pulse Time Setting](#)

## Selecting the Tray Type

### ❖ To select the tray type

1. On the Job Queue screen, select Menu by pressing the F1 key (Figure 15).

**Figure 15.** Autosampler control terminal with the Job Queue screen displayed



2. From the Menu screen, open the Tray screen as follows:
  - a. Rotate the outer knob to move the cursor bar (highlight) to **Setup**, and then press the inner knob (button).
  - b. Rotate the outer knob to move the cursor bar to **Objects**, and then press the inner knob (button).
  - c. Rotate the outer knob to move the cursor bar to **Trays**, and then press the inner knob (button).

The Trays screen lists the available trays.

3. Select a tray as follows:
  - a. Rotate the outer knob to move the cursor bar to the tray that you want, and then press the inner knob (button).

The screen now shows the Tray Type.

- b. Make the Tray Type active by pressing the inner knob (button) again.

The Tray Type is highlighted.

- c. Rotate the outer knob to select the tray type. Select **NONE** if no tray is present (rotate the outer knob until NONE appears as the Tray Type).
  - d. Press the inner knob (button) to make the selection.
4. Repeat [step 3](#) for all the trays to be used.
  5. Return to the Job Queue screen by pushing the F4 key (Home) on the control terminal.


## Changing the Pulse Time Setting

To control LC systems with a Thermo Scientific LC pump and an autosampler from the Xcalibur data system, you must change the default pulse time setting for the autosampler from 2 to 4 seconds. The pulse time is the length of time that the autosampler sends an output signal after it makes an injection. The LC pump receives this output signal through the contact closure cable.

**IMPORTANT** If your LC system contains a Thermo Scientific LC pump and an autosampler, you must change the autosampler's default pulse time setting. If you leave the default pulse time set to 2 seconds, the status of the LC pump remains at waiting for contact closure after you make an injection.

Change the pulse time setting to 4.0 seconds from the Injected screen ([Figure 16](#)).

**Figure 16.** Menu path to Injected screen



Injected	
Destination	SW-Out1
Pulse Time	4.0 s
Oms...9.9s	
Set	Home
Deflt	

### ❖ To change the pulse time setting

1. Open the Injected screen as follows:
  - a. On the Job Queue screen, press the F1 key to select **Menu** (Figure 15 on page 20).
  - b. Rotate the outer knob to move the cursor bar (highlight) to **Setup**, and then press the inner knob (button).
  - c. Rotate the outer knob to move the cursor bar (highlight) to **Objects**, and then press the inner knob (button).
  - d. Rotate the outer knob to move the cursor bar (highlight) to **Out Signals**, and then press the inner knob (button).
  - e. Rotate the outer knob to move the cursor bar (highlight) to **Injected**, and then press the inner knob (button) to display the Injected screen (Figure 16 on page 21).
2. Rotate the outer knob to move the cursor bar (highlight) to **Pulse Time**, and then press the inner knob (button) to highlight the time cell.
3. Rotate the outer knob to change the pulse time to **4.0 seconds**, and then press the inner knob (button) to set this selection as the default.
4. Press the F4 key (Home) to return to the Job Queue screen.

## Specifying the Instrument Method Parameters

To inject a sample set automatically, you must create an instrument method that contains the chromatographic conditions and mass spectrometer data acquisition settings, and then create an acquisition sequence that specifies the instrument method and vial location for each run.

This section describes how to specify the basic instrument method parameters for the autosampler. For information on how to create acquisition sequences, refer to the Xcalibur Sequence Setup view Help topic.

## Creating an Instrument Method

This section provides information about how to create an instrument method.

### ❖ To specify the instrument method parameters for the autosampler



1. Choose **Start > All Programs > Thermo Xcalibur > Xcalibur**, and then click the **Instrument Setup** icon to open the Instrument Setup window.

**Note** The Instrument Setup window might take several seconds to open. If you receive a server failure error, follow the instructions on [page xi](#).

2. Click the autosampler icon to open its Instrument Setup window.
3. Click the **Method Setup** tab to open the Method Setup page ([Figure 17](#)).

**Figure 17.** Method Setup page for the autosampler (before you select a template)

The screenshot shows the "Method Setup" window with two tabs: "Method Setup" (selected) and "Method Summary". The "Method Setup" tab contains several sections: "Template Selection" with a dropdown menu and a "Browse..." button; "Template Description" with a scrollable text area; "Syringe" with a dropdown menu; "Macro Sequence" with a scrollable text area; "Recommended Injection Volume" with a "Volume (µl)" field set to "1.000"; and "Look Ahead Injections" with a "Delay Time (mins)" field set to "0.00". On the right side, there is a "Macro Variables" section with a scrollable text area. At the bottom right, there are "Help" and "Default All" buttons.

4. Under Template Selection, click **Browse** to open the Browse for Folder dialog box, and then browse to the following folder that contains the autosampler methods (.pme extension):

*drive:\Thermo\Instruments\LC Devices\AccelaOpenAS\PAL\Methods*

5. In the Template list, select a template that represents the task that you want to do.

The autosampler comes with these templates:

- DLW Priming\_Rev01.pme
- Fast Injection Accela Open.pme
- LC-Inj DLW Fast\_Rev03.pme
- LC-Inj DLW Standard\_Rev03.pme
- Multi Step LC Injection.pme
- Priming Accela Open.pme
- Single Step LC Injection.pme
- Standard Injection Accela Open.pme
- Startup-Shutdown.pme

The Accela Open Autosampler with the Dynamic Load and Wash (DLW) option uses the following three templates in instrument setup:

- Fast Injection Accela Open.pme
- Standard Injection Accela Open.pme
- Priming Accela Open.pme autosampler (comprised of three templates):
  - Multi Step LC Injection.pme
  - Single Step LC Injection.pme
  - Startup-Shutdown.pme

If you ordered the DLW option with the autosampler, you can also select the DLW templates.




**CAUTION** For the Accela Open AS with DLW, do not use the Single Step LC Injection.pme and Multi Step LC Injection.pme templates. These are settings for the Accela Open AS without DLW and the Thermo PAL autosampler. Using these templates for the Accela Open AS with DLW could cause damage to the instrument and cause the injection to not perform any washes at all.

**IMPORTANT** If you select a DLW template, make sure the Rear Volume + Injection Volume + Front Volume + 2x Airgap is less than the syringe size. Otherwise, you receive Error 86 (Syringe Volume Out of Range).

In addition to these templates, you can create your own custom templates (see [“Using the Template Editor to Create Custom Templates”](#) on page 31). When you store a custom template in the same folder as the standard templates, the template appears in the Template list.

Figure 18 shows the macro sequence and the macro variables for the Single Step LC Injection template.

**Figure 18.** Method Setup page on the Instrument Setup window with the LC-Inj template



When you select a template, the macros used by the template appear in the Macro Sequence area. When you click on a macro listed in the Macro Sequence area, the variables used by the Macro appear in the Variables area.

- Under Variables, keep the parameters set to the defaults, or type the appropriate values for your application.
- From the Syringe list, select the syringe size that is installed in the autosampler.

The default variable settings for the Single Step LC Injection template and the default setting of 10 µL in the Recommended Injection Volume box are suitable for a 100 µL syringe. Changing the syringe size can change the allowable ranges for the recommended injection volume and the variables specified in the LC-Inj macro.

**IMPORTANT** Make sure that the syringe selection matches the size of the installed syringe.

- If the autosampler is set up to recognize the installed syringe size, an error message appears when you download methods that specify a different syringe size.
- If the autosampler is not set up to recognize the installed syringe and the specified syringe size does not match the actual syringe size, the autosampler does not inject the specified injection volume.

8. Under Recommended Injection Volume, in the Volume ( $\mu\text{L}$ ) box, type the volume of sample that you want the autosampler to inject.

The allowable injection volume range is based on the syringe size.

**IMPORTANT** For partial loop injections (variable volume), make sure that the recommended injection volume is less than half the sample loop size.

9. Under Look Ahead Injections, if the Delay Time (mins) box is available, type the amount of time that you want the autosampler to wait after it completes the current injection before it starts the next injection.

During a look ahead injection, the autosampler aspirates the sample for the next injection and waits until the end of the current run to inject the sample into the valve. When you add a delay time, the autosampler does not begin the next injection cycle immediately after completing the current injection cycle. If you enter a delay time longer than the method run time, the autosampler begins the next injection at the end of the current run. This means that adding a long delay time does not add additional run time to the method run time, but it does cancel the effect of using the Look Ahead Injections feature.

**Note** The Look Ahead Injections feature is available if you selected the Enable check box in the Look Ahead Injections area of the Pal Configuration dialog box.

10. Specify the acquisition parameters for the other devices of the LC or LC/MS instrument.

11. Save the method:

- a. From the Instrument Setup window, choose **File > Save As**.
- b. Select a folder in which to store the method, and then type a name for the method in the File Name box.
- c. Click **Save**.

The File Summary Information dialog box appears.

- d. (Optional) In the Comment box, type additional information about the method.
- e. Click **OK**.

The Xcalibur data system (or other Thermo Scientific data acquisition applications such as LCquan™) stores the method as a METH file in the specified file location.

## Instrument Setup Parameters

Use the Method Setup page (Figure 18 on page 25) to specify instrument method parameters for the autosampler.

**Table 5.** Method Setup page parameters (Sheet 1 of 2)

Parameter	Description
<b>Template Selection</b>	
Template	Shows the list of default method templates.  After you select a template, the Template list changes to a box that shows the path to the selected method template file.
Method Folder	Shows the path to the default template folder before you select a default method template.  After you select a method template (in the Template list), the Method Folder box disappears.
Browse	Opens the Browse for Folder dialog box to locate the directory that contains the method templates.  After you select a method template (in the Template list), the Browse button changes to the New Template button.
New Template	Select this button to bring back the Template list, Method Folder box, and Browse button.  After you select a method template (in the Template list), this button appears.
<b>Template Description</b>	
Template Description	Shows a description of the selected template, which appears in the Template box.
<b>Syringe</b>	
Syringe	Shows a list of the allowed syringe volumes.  Select the appropriate syringe volume.
<b>Macro Sequence</b>	
Macro Sequence	Shows the sequence of macros that make up the method displayed in the Template box.  Select a macro name to be displayed in the Macro area; its variables also appear in the Variables area.

**Table 5.** Method Setup page parameters (Sheet 2 of 2)

Parameter	Description
<b>Recommended Injection Volume</b>	
Volume (µL)	<p>When you select an instrument method for use to generate a new sequence (on the Xcalibur Home Page), the injection volume is set to what is entered in this box.</p> <p>For the LCQuan application: When “From AS” (From Auto Sampler) is specified in the LCQuan injection volume sequence cell, the autosampler uses the volume in this box when you submit the LCQuan sequence for acquisition.</p>
<b>Look Ahead Injections</b>	
Delay Time (mins)	<p>Enter the Delay Time (in minutes) between injections when using the Look Ahead feature.</p> <p>You enable the Look Ahead Injections option from the Pal Configuration dialog box (Figure 14 on page 17). When you enable this option, the next Look Ahead Injection is delayed the amount of time entered in this box. The delay time countdown begins at the completion of any post-injection steps, which prevents a sample from being held in the injection syringe for the duration of the current run. If the entered Delay Time is too long, or the current run ends prematurely, the entered Delay Time is canceled and the next injection commences in the normal manner.</p>
<b>Macro</b>	
—	<p>Shows the name of the selected macro, and the sequence number (<i>x</i>) and total number of macros (<i>y</i>) (displayed as “name”, <i>x</i> of <i>y</i>).</p> <p>This area also lists the variables associated with the macro chosen in the Macro Sequence pane.</p>
Variables	<p>Shows the variable parameters used by the selected macro.</p> <ul style="list-style-type: none"> <li>• Variables with specific allowed values (discrete variables) appear in lists. You may choose the allowed values of the discrete variables from those presented in the list.</li> <li>• Continuous variables are displayed in boxes. The allowed values for continuous variables appear when you hold the cursor over the variable’s box. You can type the value for the continuous variable in the box.</li> </ul>
Default All	This button changes the variables, used in the selected Macro Sequence, to their default values.

## Viewing the Method Summary

Use the Method Summary page to view the details of the Macro Sequence that makes up an autosampler Method or Template.

The contents of the Method Summary/Template Summary page is read-only.

❖ **To view a summary of the autosampler portion of an instrument method**


1. Open the Instrument Setup view for the autosampler (see [“Specifying the Instrument Method Parameters”](#) on page 22).
2. Choose **File > Open**, and then select the instrument method.
3. Click the **Method Summary** tab.

The Method Summary page appears with a summary of the autosampler parameters for the active instrument method ([Figure 19](#)).

## 2 Getting Started

Viewing the Method Summary

**Figure 19.** Method Summary page in the Instrument Setup window



## Creating Custom Templates and Macros

This chapter describes how to create custom templates and macros. It also includes instructions on how to define variables for macro creation. The Accela Open Autosampler provides three standard templates.

### Contents

- [Using the Template Editor to Create Custom Templates](#)
- [Testing a Custom Template](#)
- [Viewing the Template Summary](#)
- [Standard Macros](#)
- [Using the Macro Editor to Create Custom Macros](#)
- [Defining Variables](#)

## Using the Template Editor to Create Custom Templates

Use the Template Editor page of the Pal Template\Macro Editor window to make a new template (a sequence of macros). Each macro in the sequence is made up of variables that you may also edit on this page. [Table 6](#) on [page 33](#) lists the parameters for the Template Editor page.


### ❖ To create a custom template

1. Open the Template Editor page as follows ([Figure 20](#)):
  - a. Open the Instrument Setup window for the autosampler (see [“Specifying the Instrument Method Parameters”](#) on [page 22](#)).
  - b. Choose **Accela Open AS > Template\Macro Editor**.

### 3 Creating Custom Templates and Macros

Using the Template Editor to Create Custom Templates

**Figure 20.** Template Editor page on the Pal Template/Macro Editor window



2. (Optional) In the Template Description box, type a description of your custom template.
3. In the Syringe list, select the appropriate syringe size to be specified in the template.
4. Add macros to the template as follows:
  - a. Under Macro, select a macro from the list.

The Insert Above and Insert Below buttons become available.
  - b. Click **Insert Above** or **Insert Below**.

The selected macro appears in the Macro Sequence box.
5. Delete macros from the Macro Sequence box by selecting the macro and clicking **Delete**.
6. Replace a macro in the Macro Sequence box as follows:
  - a. Under Macro, select a macro from the list.
  - b. Select the macro that you want to replace in the Macro Sequence box.
  - c. Click **Replace**.
7. Save the template:
  - a. Click **Save As**.
  - b. Select a file location, and then type an appropriate name in the File Name box.
  - c. Click **Save**.

The autosampler saves the template with the .pma file extension.

**Table 6.** Template Editor page parameters (Sheet 1 of 2)

Parameter	Description
Template Description	Enter the method description in this box.
Syringe	Select the syringe volume from this list.
<b>Macro Sequence</b>	
Macro Sequence	Shows the macros selected for the method being made in the Method Template Editor.  The macros in this sequence are selected from the <a href="#">Macro</a> list. Double-clicking a macro under Macro Sequence displays it in the Macro list as described later in this table, along with the variables that make up the macro (in the Variables area).
Insert Above	Inserts a macro from the <a href="#">Macro</a> list above the macro selected under Macro Sequence.
Insert Below	Inserts a macro from the <a href="#">Macro</a> list below the macro selected under Macro Sequence.
Replace	Replaces the macro selected under Macro Sequence with the macro that is in the <a href="#">Macro</a> list.
Delete	Deletes the macro selected under Macro Sequence.
Run	Automatically executes the macros that constitute the Method (Macro Sequence).
Step	Executes the Method (Macro Sequence) in a stepwise manner.
<b>Macro</b>	
Macro	Select individual PAL macros in this list.  When you select a macro, it is highlighted in blue. The sequence number ( <i>x</i> ) of the selected macro and the total number of macros ( <i>y</i> ) in the Macro Sequence box (“Macro <i>x</i> of <i>y</i> ”) appear at the top of the Macro list. You can place the selected macro under <a href="#">Macro Sequence</a> by using the Insert Above, Insert Below, or Replace buttons as appropriate.
Default All	Changes all of the Method Template Editor variables to their default values.

**Table 6.** Template Editor page parameters (Sheet 2 of 2)

Parameter	Description
Variables	Shows the variables that make up the macro in the Macro list. <ul style="list-style-type: none"><li>• Variables with specific allowed values (discrete variables) appear in lists. You may select the allowed values of the discrete variables from those presented in the list.</li><li>• Continuous variables are displayed in boxes. The allowed values for continuous variables appear when you hold the cursor over the variable's box. You can type the value for the continuous variable in the box.</li></ul>
<b>Buttons</b>	
Save Template	Saves the Method Template under the current file name (overwrites the current file).
Save As	Saves the Method Template under a new file name or in a new location that you type in the Save As dialog box.

## Testing a Custom Template

This section provides information about testing a custom template. Use the Run Auto Sampler Method dialog box to enter the Methods information required for the autosampler to carry out autosampling. [Table 7 on page 36](#) lists the parameters for the Run Auto Sampler Method dialog box.

### ❖ To test your new template step-by-step


1. Open the Template Editor dialog box (see [“To create a custom template” on page 31](#)).
2. Open the template to review by choosing **File > Open Method Template**.

The default path to open a template is as follows:

*drive:\Thermo\Instruments\LC Devices\AccelaOpenAS\PAL\Methods*

3. Under Macro Sequence, click **Step** to open the Run Auto Sampler Method in Step Mode dialog box (Figure 21).


**Figure 21.** Run Auto Sampler Method In Step Mode dialog box



4. Make the appropriate entries and selections as follows:
  - Under Vial Index, type an appropriate vial location for the template you are testing.
  - Under Trays, select the appropriate trays.
  - In the Injection Per Vial box, type the number of injections that you want to make from each vial.
  - Under Volume, type the volume that you want to inject.
5. Click **Run**.

The autosampler performs the actions specified in the template, step-by-step. At the completion of each step, the autosampler pauses the injection sequence and prompts you with the Pal Method Editor dialog box (Figure 22).

**Figure 22.** Pal Method Editor dialog box



6. Do one of the following:
  - To continue stepping through the template, click **Step**.
  - To run the remaining portion of the template, click **Run**.
  - To stop the run, click **Stop Run**.


### 3 Creating Custom Templates and Macros

#### Testing a Custom Template

❖ **To run the complete sequence of macros in the template without pausing**

1. Open the Template Editor dialog box (see “To create a custom template” on page 31).
2. Open the custom template by choosing **File > Template > Open Method Template**.
3. Under Macro Sequence, click **Run** to open the Run Auto Sampler Method dialog box (Figure 23).

**Figure 23.** Run Auto Sampler Method dialog box



4. Make the appropriate entries and selections (see Table 7).
5. Click **Run**.

The autosampler executes the macros listed in the template.

**Table 7.** Run Auto Sampler Method dialog box parameters (Sheet 1 of 2)

Parameter	Description
<b>Vial Index</b>	
1	Type the vial index to be used with the method (SL.index).
2	Type the vial index to be used with the method (SL.index2).
<b>Trays</b>	
1	Select the tray to be used with the method (SL.tray).
2	Select the second tray to be used with the method (SL.tray2).
<b>Injectors</b>	
1	Select the injector to be used with the method (SL.injector).
2	Select the second injector to be used with the method (SL.injector2).
<b>Volume</b>	
1	Type the injection volume to be used with the method (SL.volume).

**Table 7.** Run Auto Sampler Method dialog box parameters (Sheet 2 of 2)

Parameter	Description
2	Type the second injection volume to be used with the method (SL.volume2).
<b>Injections Per Vial</b>	
—	Type the number of replicate injections to be made from each vial.
<b>Buttons</b>	
Run	Run the selected method.
Cancel	Cancel any changes and close the dialog box.

## Viewing the Template Summary

### ❖ To view a summary of a template

1. Open the Template Editor dialog box (see “[To create a custom template](#)” on [page 31](#)).
2. Open the template that you want to review by choosing **File > Open Method Template**.

The default path to open a template is as follows:

*drive:\Thermo\Instruments\LC Devices\AccelaOpenAS\PAL\Methods*


3. Click the **Template Summary** tab.

The Template Summary page appears with a summary of the active template ([Figure 24](#)).

### 3 Creating Custom Templates and Macros

Viewing the Template Summary

**Figure 24.** Template Summary page




## Standard Macros

Macros are the building blocks that you use to create templates. [Figure 25](#) shows the macros that are provided with a standard installation. The default path to the macros is as follows:

*drive:\Thermo\Instruments\LC Devices\AccelaOpenAS\PAL\Macros*

**Figure 25.** Standard macros




## Using the Macro Editor to Create Custom Macros

Use the Macro Editor page of the Pal Template\Macro Editor window when you want to create custom macros. [Table 8](#) on [page 41](#) lists the parameters for the Macro Editor page.

### ❖ To create a custom macro

1. Open the Macro Editor page as follows ([Figure 26](#)):
  - a. Open the Instrument Setup window for the autosampler (see [“Specifying the Instrument Method Parameters”](#) on [page 22](#)).
  - b. Choose **Accela Open AS > Template\Macro Editor**.
  - c. Click the **Macro Editor** tab.

**Figure 26.** Macro Editor page of the Template\Macro Editor window



2. Under Macro, select a macro from the list.

The description of the selected macro appears under Macro Description and a list of macro variables appears under Macro Variables.

3. To add a variable to the macro, do the following:

- Click **Add** to open the Variable Definition dialog box (Figure 27 on page 44).
- Define the variable (see “Defining Variables” on page 44).

4. Save the macro as follows:
  - a. Click **Save As**.
  - b. Select a file location, and then type an appropriate name in the File Name box.
  - c. Click **Save**.

The autosampler saves the macro with the .pma file extension.

**Table 8.** Macro Editor page parameters (Sheet 1 of 3)

Parameter	Description
Macro	This list displays existing autosampler macros. Select the macro that you want to edit.
Macro Description	Displays a description of the selected macro.
<b>Macro Variables</b>	
—	Displays the variables used by the selected macro.
Add	<p>Adds a new variable to the macro variables.</p> <p>For details, see “<a href="#">Defining Variables</a>” on <a href="#">page 44</a>. You define the new variable in the Variable Definition dialog box that opens when you click Add. In this dialog box, you can select the A Numeric Value or An Object option.</p> <ul style="list-style-type: none"> <li>• If you select a numeric variable type, you see a pane for entering the new variable name, a list for the variable units, and panes for entering the default value for the variable along with the lower and upper limits for the variable.</li> <li>• If you select an object variable type, you see a check list of objects that the variable might represent.</li> </ul>
Edit	<p>Edits an existing macro variable.</p> <p>You edit a variable in the Edit Parameter dialog box that opens when you click Edit. In this dialog box, you can select the A Numeric Value or An Object option.</p> <ul style="list-style-type: none"> <li>• If you select a numeric variable type, you see a pane for entering a new variable name, a list for the variable units, and panes for entering the default value for the variable along with the lower and upper limits for the variable.</li> <li>• If you select an object variable type, you see a check list of objects that the variable might represent.</li> </ul>
Delete	Deletes the selected Macro Variable.

**Table 8.** Macro Editor page parameters (Sheet 2 of 3)

Parameter	Description
<b>Command Sequence</b>	
—	Lists the sequence of commands that make up the macro selected in the Macro list.
	When you select a command, the command parameters and description appear in the Command area. The command parameters appear in both boxes (continuous variables) and lists (discrete variables).
Insert Above	The command shown in the Command list is inserted above the command selected (highlighted in blue) under Command Sequence.
Insert Below	The command shown in the Command list is inserted below the command selected (highlighted in blue) under Command Sequence.
Replace	The command selected (highlighted in blue) under Command Sequence is replaced by the command shown in the Command list.
Delete	Deletes the command selected (highlighted in blue) under Command Sequence.
<b>Command</b>	
—	Lists the individual autosampler commands. These individual commands make up the command sequence of an autosampler macro. Select the command by moving the cursor over the command (no click is needed), highlighting each command in blue. Clicking the command displays it in this list (highlighted in blue). Simultaneously, the command parameters appear in the Parameters area below the Command list.
	The pane below the Parameters area describes the function of the selected command.

**Table 8.** Macro Editor page parameters (Sheet 3 of 3)

<b>Parameter</b>	<b>Description</b>
Parameters	Shows the command parameters used by the selected autosampler command. <ul style="list-style-type: none"> <li>Parameters with specific allowed values (discrete parameters) appear in lists. You may choose the allowed values of the discrete parameters from those presented in the list.</li> <li>Continuous parameters are displayed in boxes. The allowed values for continuous parameters appear when you hold the cursor over the respective boxes for the parameter. You can type the value for the continuous parameter in its boxes and view other allowed continuous parameter values in a shortcut menu. To activate this menu, right-click the continuous parameter box of interest. The parameter you select from the shortcut menu appears in the box.</li> </ul>
Command Description	Shows a description of the selected command in the bottom pane.
<b>Buttons</b>	
Save Macro	Saves the macro (command sequence) under the current file name (overwrites the current file).
Save As	Saves the macro (command sequence) under a new file name, in a new location, or with both a new file name and location that you type in the Save As dialog box.


## Defining Variables

Variables are the building blocks that you use to create macros. Use the Variable Definition dialog box to specify the details of a numerical or object macro variable that you want to add or edit in the Macro Variables box. Table 9 on page 45 lists the parameters for the Variable Definition dialog box.

#### ❖ To create custom variables


1. Open the Macro Editor page (“Using the Macro Editor to Create Custom Macros” on page 39).
2. Click **Add** to open the Variable Definition dialog box (Figure 27).

**Figure 27.** Variable Definition dialog box



3. Select a variable type:
  - To create a variable that is a numeric value, select the **A Numeric Value** option, and then go to [step 4](#).
  - To create an object variable, select the **An Object** option. A list of objects appears. Go to [step 5](#).
4. Complete the definition of a numeric variable as follows:
  - a. In the Variable Name box, type an appropriate name.
  - b. In the Variable Units list, select a unit of measure.
  - c. In the Default Value box, type a default value for the variable.
  - d. In the Lower Limit box, type a lower limit for the variable.
  - e. In the Upper Limit box, type an upper limit for the variable.
  - f. Click **OK** to save the variable and close the Variable Definition dialog box.
5. Complete the definition of an object variable as follows:
  - a. In the Variable Name box, type an appropriate name.
  - b. Select the appropriate objects from the list (Figure 28).
  - c. Click **OK** to save the variable and close the Variable Definition dialog box.

**Figure 28.** Object variable selections



**Table 9.** Variable Definition dialog box parameters (Sheet 1 of 2)

Parameter	Description
<b>Variable Type</b>	
A Numeric Value	Select this option if the parameter is a numeric value.  If you are editing an existing numerical parameter, this option is preselected.
An Object	Select this option if the parameter is an object.  If you are editing an existing object parameter, this option is preselected.
Variable Name	If you are editing an existing variable, this box shows the name of the variable. If you are adding a variable, type its name.
OK	Accepts the new Variable Definition or accepts the edits to an existing Variable Definition.
Cancel	Deletes changes made in the Variable Definition dialog box.
<b>Variable Type–Numeric</b>	
Variable Units	If you are editing an existing variable, this list shows the units for the selected variable. If you are adding a variable, you may enter units appropriate for that variable. When you click the list arrow, you can select the appropriate units from the options in the list.
Default Value	The nominal or center value of the variable.
Lower Limit	The lowest acceptable value of the variable.
Upper Limit	The highest acceptable value of the variable.

### 3 Creating Custom Templates and Macros

#### Defining Variables

**Table 9.** Variable Definition dialog box parameters (Sheet 2 of 2)

<b>Parameter</b>	<b>Description</b>
<b>Variable Type–Object</b>	
Please Check What Objects This Variable May Be	Select the check boxes for one or more objects that the variable describes.

# Using the Direct Controls to Operate the Autosampler

This chapter describes how to operate the Accela Open Autosampler from the Direct Control dialog box that you access from the Xcalibur Instrument Setup window or from the Inlet Direct Control dialog box that you access from the mass spectrometer Tune application. Use the direct controls when you want to control the autosampler before starting a run acquisition.

## Contents

- [Opening the Direct Control Dialog Box](#)
- [Making Single Injections from the Tune Window](#)
- [Changing the Tray Configuration](#)
- [Specifying the Location of the Direct Control Methods](#)
- [Wrapping Direct Control Injections with Additional Macros](#)

## Opening the Direct Control Dialog Box

The direct control pages for the autosampler are available from the Instrument Setup window and the Tune window.


### ❖ To open the Direct Control dialog box from the Instrument Setup window

1. In the Xcalibur Instrument Setup window, click the autosampler icon.
2. Choose **Accela Open AS > Direct Control** to open the Direct Control dialog box (Figure 29).

## 4 Using the Direct Controls to Operate the Autosampler

### Making Single Injections from the Tune Window

**Figure 29.** Direct Control dialog box



❖ **To open the Inlet Direct Control dialog box from the Tune window**

1. In the Tune window, choose **Setup > Inlet Direct Control**.
2. Click the **Accela Open AS** tab to open the Inlet Direct Control dialog box.

The Inlet Direct Control dialog box has the same functionality as the Direct Control dialog box shown in [Figure 29](#).

## Making Single Injections from the Tune Window

You can make single injections from the Tune window by using the Methods page of the Inlet Direct Control dialog box. [Table 10](#) lists the parameters for the Methods page of the Inlet Direct Control dialog box.

❖ **To make single injections from the Tune window**

1. Open the Tune window for your mass spectrometer.
2. Open the Inlet Direct Control dialog box (see [“Opening the Direct Control Dialog Box”](#) on [page 47](#)).
3. Specify the location of the methods that you want to use by using the Direct Control Configuration page (see [“Specifying the Location of the Direct Control Methods”](#) on [page 52](#)).

4. Specify the tray types that are loaded in the autosampler tray holders (see “[Changing the Tray Configuration](#)” on page 50).

5. Click the **Methods** tab.

The Methods page appears. The Syringe area lists the size of the configured syringe. The Device Status area provides the status of the autosampler.

6. In the Method list, select a method.

A description of the method appears in the Description box. The method list depends on the file locations specified on the Direct Control Configuration page.

**Note** The visible DLW methods in the Method list are only available if you installed the DLW option that comes with the autosampler. See [Chapter 5, “Using Dynamic Load and Wash \(DLW\).”](#)

7. Make the following selections and entries:

- Under Vial Index, type the location of the sample vial.
- Under Trays, select the tray where the sample vial is located.
- Under Volume, type the volume that you want the autosampler to inject.

8. In the Inlet Direct Control dialog box, click the tab for the configured LC pump and start the solvent flow from the LC pump.

9. Return to the autosampler Methods page, and then click **Run**.

**Table 10.** Methods page parameters (Sheet 1 of 2)

Parameter	Description
Method	Use this list to select a method file. The visible DLW methods in the Method list are only available if you installed the DLW option that comes with the autosampler.
Description	Contains a description of the selected method.
Syringe	Lists the size of the configured syringe.
Device Status	Displays the status of the autosampler.
<b>Vial Index</b>	
1	Type the vial index to be used with the method (SL.index).
2	Type the second vial index to be used with the method (SL.index2).

**Table 10.** Methods page parameters (Sheet 2 of 2)

Parameter	Description
<b>Trays</b>	
1	Select the tray to be used with the method (SL.tray).
2	Select the second tray to be used with the method (SL.tray2).
<b>Volume</b>	
1	Type the injection volume to be used with the method (SL.volume).
2	Type the second injection volume to be used with the method (SL.volume2).
<b>Injectors</b>	
1	Select the injector to be used with the method (SL.injector).
2	Select the second injector to be used with the method (SL.injector2).
Injections Per Vial	Type the number of replicate injections to be made from each vial.
Run	Runs the selected method.

## Changing the Tray Configuration

Use the Tray Configuration page of the Direct Control dialog box to specify the tray types to be used for direct control injections. The autosampler uses this tray information instead of the tray information specified in the autosampler Control Terminal.

Table 11 lists the parameters for the Tray Configuration page of the Direct Control dialog box. For information about the sample trays, see [Appendix A, “Sample Trays.”](#)


**IMPORTANT** The autosampler does not store the tray configuration that you specify on the Tray Configuration page of the Direct Control dialog box. When you close the dialog box, the tray configuration resets to the autosampler Control Terminal settings.

### ❖ To temporarily change the tray configuration

1. Open the Direct Control dialog box (see [“Opening the Direct Control Dialog Box”](#) on [page 47](#)).
2. Click the **Tray Configuration** tab to open the Tray Configuration page ([Figure 30](#)).

For information about the sample trays, see [Appendix A, “Sample Trays.”](#)

**Figure 30.** Tray Configuration page (VT54 example)



3. In the Tray list, select the tray that you want to modify.
4. In the Type list, select the tray type.

**Table 11.** Tray Configuration page parameters

Parameter	Description
Tray	Select the tray location in the three-drawer (deep drawer) configuration. For information about the sample trays, see <a href="#">Appendix A, “Sample Trays.”</a>  Range: Tray01–Tray12
Type	Select the tray type.
Tray Holder	The allowed values are as follows: <ul style="list-style-type: none"> <li>• Stack1, for Tray01–Tra06</li> <li>• CStack1, for Tray07–Tray12</li> </ul> Tray holders can hold one or more trays. A Stack is a tray holder that is designed to hold micro plates. The autosampler firmware determines the type and number.
Vials in X Direction, Vials in Y Direction	Determined by choice of tray type.

## Specifying the Location of the Direct Control Methods

Use the Direct Control Configuration page of the Direct Control dialog box to specify the location of the methods that you can open on the Methods page. Table 12 lists the parameters for the Direct Control Configuration page.


### ❖ To specify the method folders

1. Open the Direct Control dialog box (see “Opening the Direct Control Dialog Box” on page 47).
2. Click the **Direct Control Configuration** tab to open the Direct Control Configuration page (Figure 31).

The default path to the methods that appear in the Method list on the Method page is as follows:

*drive:\Thermo\Instruments\LC Devices\AccelaOpenAS\PAL\DirectControlMethods*

**Figure 31.** Direct Control Configuration page of the Direct Control dialog box



3. To add additional file locations, type the folder location in the Path box or click **Browse** to find the appropriate folder, and then click **Add**.

The new folder location appears in the box below the row of buttons.

4. To delete a file location, select the folder path, and then click **Delete**.
5. To reset the file location to the default folder, click **Reset**.

**Table 12.** Direct Control Configuration page parameters

Parameter	Description
<b>Wrap Each Direct Control Method Run with the Following Macros</b>	
Sequence Entry Exit	Select this check box to wrap methods with the sequence list entry and exit macros (SLEntry and SLExit).
Method Entry Exit	Select this check box to wrap methods with the method entry and exit macros (Method Entry and Method Exit).
<b>Note</b> For more information, see <a href="#">Wrapping Direct Control Injections with Additional Macros</a> .	
<b>Method Search Paths</b>	
Path	Type the extra path that you want to add for the method search routine.
Browse	Use this button to search for the path to a file folder.
Add	Use this button to add a file folder for the method search routine.
Delete	Use this button to delete a file folder from the method search routine.
Reset	Use this button to remove all of the user-added file folders.
[text box]	Displays the file folders (directories) where the autosampler searches for methods.
Firmware Version	Specifies the firmware version downloaded from the autosampler.

## Wrapping Direct Control Injections with Additional Macros

Use the Direct Control Configuration page of the Direct Control dialog box to specify whether the autosampler is to run the SL Entry and SL Exit, the Method Entry and Method Exit, or all of these macros before and after an injection.

❖ **To specify the use of the SL Entry and SL Exit macros**

Select the **Sequence Entry Exit** check box.

❖ **To specify the use of the Method Entry and Method Exit macros**

Select the **Method Entry Exit** check box.



## Using Dynamic Load and Wash (DLW)

This chapter describes how to use the Dynamic Load and Wash (DLW) option, which includes how to install the Cycle Composer Macros or ICC Cycles, how to operate the DLW, and a step-by-step illustration of the DLW cycles.

### Contents

- [Overview](#)
- [Cycle Composer Macros or ICC Cycles](#)
- [Operating Dynamic Load and Wash \(DLW\)](#)
- [DLW Cycle Step-by-Step](#)

## Overview


This section contains an overview of the Dynamic Load and Wash (DLW) option.

The DLW option ([Figure 32](#)) represents a new wash station concept that combines an injection cycle with wash steps. The linked combination of the two steps, which are usually separate, minimizes cycle time and carryover.

There are two characterizing features for the DLW option:

- The sample solution never contacts the syringe itself; it is held between the holding loop.
- Wash solvents are pumped from back to front into the DLW system to intensely flush all critical parts that are in contact with the sample.

**Figure 32.** Dynamic Load and Wash (DLW)



The DLW option consists of two self-priming micro pumps (mounted on a dedicated bracket), and the wetted parts are Ryton™ PPS and Kalrez™ (FFPM). The pump IN ports connect to the wash solvent bottles and the OUT ports connect to the DLW manifold, which is part of the assembly of the dedicated DLW syringe holder. A holding loop separates the syringe from the DLW actuator to prevent sample from contacting these parts.

The syringe and holding loop are preloaded with wash solvent #1 at the start. The sample is picked up and remains separated from wash solvent #1 by an air gap. After loading the loop and injection, wash solvent #1 is pushed into the system, followed directly by wash solvent #2 to flush the critical valve paths.

The DLW syringe assembly moves to the wash station for further cleaning steps and for preparing the syringe and holding loop for the next cycle.

For further details, see “DLW Cycle Step-by-Step” on [page 68](#).

## Cycle Composer Macros or ICC Cycles

You can only operate the autosampler DLW option with PAL control software, such as the ThermoPAL driver for the Xcalibur data system or Cycle Editor for ICC interpretation in another data handling system software.

For software control, the autosampler provides three macros or three cycles. The first macro covers the initial and daily priming of the solvent lines and covers a complete injection cycle. The second macro is used for Standard speed with optimized washing possibilities. The third is used for a Fast cycle for optimized throughput and less focus on carryover. See [Table 13](#).

**Table 13.** The DLW macro definitions

Macro name	Macro description
Priming Accela Open	For initial and daily routine priming of the solvent lines and DLW manifold. The Clean Time for both wash solvents is a variable for the user to define the intensity of washing.
Standard Injection Accela Open	Standard injection cycle using all possibilities of the DLW option. The injection valve inlet port and the needle are washed with both wash solvents (inside and out). You can add an extra Stator Wash for intensive washing of the injection valve (valve toggle).
Fast Injection Accela Open	Tuned for speed and high-throughput application. It differs from the Standard DLW macro in that some steps are left out to shorten the cycle time.

For detailed descriptions of all macros (or cycles), see [Table 14](#) on [page 59](#) and [Table 15](#) on [page 60](#).

**Note** The provided macros are written for standard injection valve drives, which are controlled and activated through the AUX interface.

# Installing the Cycle Composer Macros or ICC Cycles

The autosampler ships with a DVD containing various cycles for the DLW option. Macros for the Xcalibur data system are installed by the LC Devices installer. Copy these macros to the driver method folder or the corresponding folder for application within the integrated system.

### ❖ To copy cycles to your Cycle Composer

1. Navigate in Windows Explorer to the ThermoPAL driver folder. The usual folder location is as follows:

*drive:\Thermo\Instruments\LC Devices\AccelaOpenAS\PAL\*

2. If you want to add the DLW option macros to an existing method folder, copy the macro (.pma) and the method (.pme) files from the DLW Option folder on the DVD to the method folder.

If the PAL system is integrated in a data system software that controls the autosampler by using the Cycle Editor for PAL ICC interpretation (for example, Analyst™, ChemStation™, Empower™, EZChrom, MassLynx™, or Xcalibur), an ICC Cycle is used and not the Cycle Composer Macro. The cycle file extension is .cyx.

**Tip** You can convert a Cycle Composer macro to a cycle (extension .cyx) by using the Cycle Editor software. Conversion is available starting with Cycle Editor version 1.4.0.4.

## General Considerations

You must establish the duration of the wash steps for each configuration and application. Consider factors such as the viscosity and surface tension of the individual wash solvent composition and the backpressure of the system.

Be aware that a higher backpressure builds up if the valve bore size (standard valve bore 0.25 mm) or the installed loop internal diameter is lower. Standard loop internal diameter (ID) for Thermo-defined loops with a volume of 5, 10, and 20 µL is 0.25 mm. The loop with 2 µL content volume has an ID of 0.125 mm.

Keep the tubing internal diameters of the tubing in line with the valve dimensions, loop ID, and flow rate.

## Priming the Solvent Lines, Wash1 and Wash2

Use this macro at installation to prime the entire system. Set the wash time to approximately 120 seconds for each solvent.

After installation, for best results prime the system before activating the first run. For daily preparation of the system, the wash time can be much shorter: approximately 20 seconds. The goal is an entire liquid system free of any air bubbles.

### ❖ To prime the solvent lines

1. Set the wash time to approximately 120 seconds for each solvent.
2. After installation, for best results prime the system before activating the first run. For daily preparation of the system, the wash time can be much shorter: approximately 20 seconds. The goal is an entire liquid system free of any air bubbles.

Table 14 describes the function of the Priming Accela Open AS\_Rev01 macro.

**Table 14.** Priming Accela Open AS\_Rev01 macro (Sheet 1 of 2)

Macro description	Macro variable
The injection unit moves to the DLW Wash Station, position Wash1.	Clean Time solvent 1 Eject Speed DLW Syringe
The injection unit moves to the DLW Wash Station, position Wash2.	Clean Time solvent 1 Eject Speed DLW Syringe
The DLW system is rinsed with Wash Solvent 1 in position Waste.	—
The injection valve is cleaned first with the content of the DLW syringe (Wash Solvent 1), followed by Wash Solvent 2 and finally, the last wash to prepare the system for an injection cycle, rinsed with Wash Solvent 1 again.	Needle Gap is a parameter from Rinse Inj Atom. The variable in this macro is Needle Gap Valve Clean.  The function of this parameter is to minimally raise the needle in the injection port to allow rinsing around the needle tip.
Remark: The Atom Rinse Inj is new, available starting with FW 4.1.x. The DLW actuator/solenoid is activated; the Wash Solvent (pump), the Needle Gap, and the Rinse Time are selectable.	The pressure of the spring-loaded balls in the DLW Syringe Holder assembly is released by moving approximately 3 mm up (default). This leaves a gap, between the needle tip and the valve bottom, of approximately 1 mm to enable a flush at this contact point.
A Repeat-End loop enables adding an extra rinsing step, valve toggle.	Stator Wash: Counter 0 = disable valve toggle steps Counter 1 = enable valve toggle steps

**Table 14.** Priming Accela Open AS\_Rev01 macro (Sheet 2 of 2)

Macro description	Macro variable
If the counter is set to 1, follow the described steps below. If the counter is set to 0, the macro finishes at this point.	—
The injection moves to the injection valve. The valve is switched to the Active position.	—
The valve is rinsed with Wash Solvent 2, followed by Wash Solvent 1.	Stator Wash Time Solvent 2 Stator Wash Time Solvent 1  Remark: The loop is filled with the last rinse of Wash Solvent 1. Verify the composition of Wash Solvent 1. The solvent should have a lower elution power than the solvent gradient starting conditions or sample solvent composition. This is important for partial loop filling.
The injection valve is switched back to the Standby position.	—
End of macro DLW Priming.	—

## Standard DLW Injection Cycle

Table 15 describes the function of the Stand Injection Accela Open macro.

**Table 15.** Stand Injection Accela Open macro parameters (Sheet 1 of 3)

Macro description	Macro variable
The PAL system waits first for the Sync Signal Ready before the injection cycle starts.	Remark: Sync Signal setting Start
The injection valve moves to a defined Standby position.	Inject to Standby
The Rear air segment pulls into the Holding Loop.	Airgap Volume Filling Speed
The sum of Rear-, Sample List-, and Front-Volume aspirates into the Holding Loop.	Front Volume Rear Volume (SL.volume)
The Front air segment aspirates.	Airgap Volume Filling Speed Pullup Delay

**Table 15.** Stand Injection Accela Open macro parameters (Sheet 2 of 3)

Macro description	Macro variable
The injection unit moves to the DLW Wash Station, Wash1 position.	—
The needle is inserted (dipped) for 1 second to wash the outer needle surface. No plunger movement at this step.	
The injection unit moves to the specified injection valve. The Front- and Airgap-Volume is ejected.	Inject to Front Volume Airgap Volume Injection Speed
The PAL system waits for the data system.	Wait for DS
The injection valve switches to Active position. Wait the Pre Inject Delay time.	Inject to Pre Inject Delay
The loop fills with the sample volume as specified in the Sample list.	(SL.volume) Injection Speed
The injection valve switches to Standby position, the loop content is injected. Timer Delay Stator Wash starts and sends a Start signal to the HPLC system.	Post Inject Delay Timer 1
The plunger of the DLW Syringe moves down to dispense the Rear Sample and Air Segment to Waste. The Holding Loop is still filled with Wash Solvent 1.	(Syr. Eject Speed)
The DLW Actuator/Solenoid activates to deliver Wash Solvent 2 into the Holding Loop to clean the injection valve from Port 1 to Port 2.	Wash2 Inject to Needle Gap Valve Clean Valve Clean Time Solvent 2
For this step the needle tip is lifted, releasing the sealing pressure to enable rinsing around the tip sealing point.	
The injection unit moves to the DLW Wash Station, Wash2 position.	Wash2 (Syr.Eject Speed)
The needle is rinsed inside and out with Wash Solvent 2.	Post Clean Time Solvent2

## 5 Using Dynamic Load and Wash (DLW)

Cycle Composer Macros or ICC Cycles

**Table 15.** Stand Injection Accela Open macro parameters (Sheet 3 of 3)

Macro description	Macro variable
The injection unit moves back to the injection valve. The Inlet Port and engraving to waste Port are flushed with Wash Solvent 1 to prepare the valve for the next injection.	Wash 1 Inject to Needle Gap Valve Clean Valve Clean Time Solvent 1
The injection unit moves back to the DLW Wash Station, Wash1 position to flush the syringe needle inside and out with Wash Solvent 1.	Wash 1 Post Clean Time Solvent 1
This is a preparation step for the next injection, and especially important for biofluid samples.	
Cycle end for LC-Inj DLW Standard macro.	—
An optional cleaning step is attached to the DLW Standard injection cycle: Stator Wash or valve toggle.	Stator Wash Stator Wash count: 1 = Cleaning step active Stator Wash count: 0 = Cleaning step disabled
A Repeat-End loop can be activated with the Count.	
If Stator Wash is activated, the following steps are executed.	—
The injection unit moves to the injection valve. From the last step above, the Holding Loop is filled with Wash Solvent 1.	Inject to Delay Stator Wash (Active)
Timer 1 waits to switch the valve (Toggle) into Active position (fill loop).	
The DLW Actuator/Solenoid activates to deliver Wash Solvent 2 to the Holding Loop and into the valve system.	Inject to Wash2 Stator Wash Time Solvent 2
The first solvent flush arriving at the valve is Wash Solvent 1 parked in the Holding Loop at the beginning, followed by Wash Solvent 2.	
Wash Solvent changes to Wash Solvent 1.	Inject to Wash 1 Stator Wash Time Solvent 1
The injection valve switches back to the Standby position.	Inject to (Standby)

## Fast DLW Injection Cycle

The Fast injection cycle differs from the Standard cycle as follows:

- The needle is not dipped in the Wash station Wash1 after sample pickup and before it moves to the injection valve.
- The wash steps after injection are reduced to Valve Clean with Wash Solvent 1 and Wash Solvent 2. The DLW needle is flushed in the DLW Wash Station with Wash Solvent 1 only.
- Stator Wash (valve toggle) is not available.

Table 16 describes the function of the Fast Injection Accela Open macro.

**Table 16.** Macro Fast Injection Accela Open (Sheet 1 of 2)

Macro description	Macro variable
The PAL system waits first for the Sync Signal Ready before the injection cycle starts.	Remark: Sync Signal setting Start
The injection valve moves to a defined position: Standby.	Inject to Standby
The Rear air segment aspirates into the Holding Loop.	Airgap Volume Filling Speed
The sum of Rear-, Sample List-, and Front-Volume aspirates into the Holding Loop.	Front Volume Rear Volume (SL.volume)
The Front air segment aspirates.	Airgap Volume Filling Speed Pullup Delay
The injection unit moves to the specified injection valve. The Front- and Airgap-Volume ejects to Waste.	Inject to Front Volume Airgap Volume Injection Speed
The PAL system waits for the data system.	Wait for DS
The injection valve switches to Active position. Wait the Pre Inject Delay time.	Inject to  Pre Injection Delay
The loop fills with the sample volume as specified in the sample list.	(SL.volume)
The injection valve switches to Standby position, and the loop content is injected.	Injection Speed Post Inject Delay

## 5 Using Dynamic Load and Wash (DLW)

### Operating Dynamic Load and Wash (DLW)

**Table 16.** Macro Fast Injection Accela Open (Sheet 2 of 2)

Macro description	Macro variable
The plunger of the DLW Syringe moves down to dispense the Rear Sample and Air Segment to Waste. The Holding Loop is still filled with Wash Solvent 1.	Injection Speed
The DLW Actuator/Solenoid activates to deliver Wash Solvent 2 into the Holding Loop to clean the injection valve from Port 1 to Port 2.	Inject to Wash2 Needle Gap Valve Clean Valve Clean Time Solvent 2
For this step the needle tip is lifted, releasing the sealing pressure to enable rinsing around the tip sealing point.	
Wash Solvent 1 follows to prepare the valve for the next injection.	Wash1 Inject to Needle Gap Valve Clean Valve Clean Time Solvent 1
The injection unit moves back to the DLW Wash station, Wash1 position to flush the syringe needle inside and out with Wash Solvent 1.	Wash1 Post Clean Time Solvent 1
This is a preparation step for the next injection, and is especially important for biofluid samples.	
Cycle ends for LC-Inj DLW Fast macro.	—

## Operating Dynamic Load and Wash (DLW)

This section describes how to operate the Dynamic Load and Wash (DLW) option.

- [Priming the Solvent Lines](#)
- [Location of Solvent and Waste Bottles](#)
- [DLW Pumps](#)

## Priming the Solvent Lines

**Tip** For trouble-free DLW operation, make sure the two solvent lines are free of air bubbles at all times. If the solvent lines are being connected for the first time or during a solvent change, you must prime the solvent lines properly until air bubbles are no longer visible. Use solvent degassing for best results.

To make the initial and daily priming efficient and controllable, the Open autosampler comes with a Cycle Composer macro, or ICC cycle.

### ❖ To prime the solvent lines

1. Load the macros and methods into the folder.
2. Start the corresponding macro for initial or daily priming.
3. Check the solvent lines and prime until air bubbles are no longer visible.
4. Press F4 for Home.

## Location of Solvent and Waste Bottles

The DLW option contains self-priming membrane pumps. You can place the solvent bottles either in the fast wash station holder or on the lab bench.

You must place the Waste bottle greater than 30 cm (11.8 in.) below the injection valve. Make sure that the waste liquid can flow into the waste bottle without restriction. Place the waste tubing above the level of the liquid. Ideally, the tube is fixed at the neck of the waste bottle.

**Tip** Use good laboratory practice to avoid contaminating the wash solvents and the wash bottles. Avoid biological growth in pure water by either replacing it regularly or adding a small percentage of organic solvents, such as methanol or acetonitrile. Certain buffer solutions can decompose at room temperature when exposed to light. Filtering the wash solvents before filling the bottle, especially if using salt buffers, is mandatory to avoid any clogging of the solvent paths.

### DLW Pumps

From the control point of view, the DLW pumps respond in the same manner as the fast wash station. Power-out signals activate the pumps. Because the electric current setting for the DLW is different, you must load the corresponding PAL Firmware Objects for the DLW wash station type.

The wetted parts in the pump are made from the following materials:

- Membrane: Kalrez (FFPM)
- Body, valves: Ryton PPS

The pumps are self priming with a maximum suction lift of a 3 m column of water.

### DLW Actuator/Solenoid

The DLW Actuator/Solenoid has the function of separating and completely shutting off the lines in the direction of the syringe (sample loading) or the wash solvent lines.

After opening the DLW Actuator/Solenoid for the wash solvent lines, you can pump the desired wash solvent into the system by activating the corresponding DLW pump.

[Figure 64](#) on [page 98](#) illustrates this functionality as part of the step-by-step cycle for a standard injection.

The wetted parts in the DLW Actuator Solenoid are made from the following materials:


- Solenoid body: PEEK™
- Seal material: FFKM (Simriz™)

**Note** PEEK exhibits excellent chemical resistance to most of the chemicals used. However, the following solvents are not recommended for use with PEEK: DMSO, THF, methylene chloride (dichloromethane), nitric acid, or sulfuric acid. For more details, refer to the compatibility tables provided by the manufacturer of PEEK material or components.


**Note** Current applied from the actuator control PCB to the actuator/solenoid activates a green LED. This activation does not indicate that the solenoid opens or closes.

[Figure 33](#) shows the DLW manifold and the actuator/solenoid connections. [Figure 34](#) shows the DLW flow diverter being inserted.

**Figure 33.** DLW manifold and Actuator/Solenoid



**Figure 34.** Inserting the DLW Flow Diverter



## DLW Cycle Step-by-Step


This section provides illustrations to demonstrate a step-by-step DLW cycle.

- [Cycle for Standard Injection](#)
- [Additional Valve Toggle Step to DLW Standard Cycle](#)
- [Cycle for Fast Injection](#)


### Cycle for Standard Injection

Figure 35 to Figure 48 illustrate a step-by-step cycle for the standard injection.

**Figure 35.** Standard: Start the cycle



**Figure 36.** Standard: Step 1 – Aspirate rear air segment



**Figure 37.** Standard: Step 2 – Get sample and aspirate rear, inject, and front volume







Figure 38. Standard: Steps 3 - 4 – Aspirate front air segment




**Figure 39.** Standard: Steps 5 - 6 – Insert the needle in wash position 1 to clean the outside




**Figure 40.** Standard: Steps 7 - 8 – Dispense front air segment and front sample volume to waste

**Figure 41.** Standard: Steps 9 - 10 – Switch valve to LOAD position and fill loop with inject volume




**Figure 42.** Standard: Step 11— Switch valve to INJECT position and start chromatographic process




## 5 Using Dynamic Load and Wash (DLW)

DLW Cycle Step-by-Step

**Figure 43.** Standard: Step 12 – Dispense rear sample volume and air segment to waste



**Figure 44.** Standard: Steps 13 - 14 – Clean valve with Wash Solvent 2



## 5 Using Dynamic Load and Wash (DLW)

DLW Cycle Step-by-Step

**Figure 45.** Standard: Steps 15 - 16 – Wash the syringe needle with Wash Solvent 2





Figure 46. Standard: Steps 17 - 18 – Clean valve with Wash Solvent 1



## 5 Using Dynamic Load and Wash (DLW)

DLW Cycle Step-by-Step

**Figure 47.** Standard: Steps 19 - 20 – Wash the syringe needle with Wash Solvent 1





Figure 48. Standard: End cycle



## Additional Valve Toggle Step to DLW Standard Cycle

This section contains information about additional steps that are necessary for a DLW Standard Cycle.

### Considerations for Additional Stator Wash Cleaning Step

The DLW Standard Cycle has the built-in option for the user to switch the injection valve at the end of the chromatographic run before equilibration of the column to the start conditions.


If the method variable “Stator Wash” is set to “1,” the extra cleaning process for the valve, with “Valve Toggle,” is part of the standard cycle.

If the method variable is deactivated (setting “0”), the DLW Standard cycle ends as shown in [Figure 48](#) on [page 81](#).

The macro (cycle) is written such that the optional valve toggle steps can be executed before re-equilibration of the column. You must synchronize the time to switch the valve with the chromatographic method by using the method variable Delay Stator Wash. The two wash solvents are timed by the method variables, Stator Wash Time Solvent 1 and Stator Wash Time Solvent 2. After these wash times have elapsed, the valve switches back to the start position.

[Figure 49](#) illustrates the recommended retention time for Stator Wash or Valve Toggle times.

**Figure 49.** Timing for Stator Wash Step



From a chromatographic viewpoint, the optional cleaning step is important to understand. Assuming that the valve stator between ports 1 and 6 (for example, in the standard Cheminert™ valve) is contaminated and cannot be cleaned during the injection process, the valve toggle brings the engraving back between the two ports. Flushing the valve with both wash solvents eliminates remaining sample material located between stator ports 1 and 6.

What points must you consider when you use the Stator Wash or Valve Toggle option?

Observe the rules if biofluid samples are injected. First, sample contact should always be with an aqueous solution to avoid protein precipitation. After washing with organic solvent (higher elution power), flush the system again with wash solvent 1.

The first toggle near the end of the chromatographic cycle provides the advantage that the sample loop is already flushed out first with the mobile phase with a solvent of high elution power (assuming gradient application).

The second valve toggle time follows immediately after finishing the second solvent flush. You cannot program a second switching time. The waiting time for the second valve toggle should be long enough so that the entire system is flushed out by both wash solvents.

Consider the entire delay volume to determine the second valve switch. The DLW internal volumes are as follows:

- Manifold, 90  $\mu$ L
- Holding Loop, 108  $\mu$ L
- Syringe Needle Gauge 22, 6.7  $\mu$ L
- Installed Injection Loop

Total delay volume: 205  $\mu$ L + Loop content volume

Do the second valve toggle (back to starting condition) before the system equilibration time has started. The Loop content is ideally a solvent of a low elution power when switched back.


If you apply isocratic chromatography, the remaining contaminants might be washed into the system and can build up higher background noise for the column, the detector, or both over a longer period of time.

## Additional Cleaning by Using the Stator Wash or Valve Toggle


Figure 50 to Figure 55 illustrate additional cleaning steps for “Stator Wash.”

### Stator Wash: End of Standard Injection Cycle


**Figure 50.** Start the additional Valve Toggle cleaning step




**Figure 51.** Stator Wash: Step 1 - 2 – Switch valve to Load Position (toggle) and clean valve with Wash Solvent 1




**Figure 52.** Stator Wash: Step 3 – Clean valve with Wash Solvent 2




**Figure 53.** Stator Wash: Step 4 – Wash Solvent 1 dispenses Wash Solvent 2



**Figure 54.** Stator Wash: Step 5 – Wash Solvent 1 cleans second valve




**Figure 55.** Stator Wash: Step 6 – Switch valve back to Inject Position (toggle)




## Cycle for Fast Injection

Figure 56 to Figure 67 illustrate a step-by-step cycle for the Fast Injection.

**Figure 56.** Fast: Cycle starts



**Figure 57.** Fast: Step 1 – Aspirate rear air segment



## 5 Using Dynamic Load and Wash (DLW)

### DLW Cycle Step-by-Step

**Figure 58.** Fast: Step 2 – Get sample and aspirate rear, inject, and front volume







Figure 59. Fast: Step 3 - 4 – Aspirate front air segment




**Figure 60.** Fast: Steps 5 - 6 – Dispense front air segment and front sample volume to waste




**Figure 61.** Fast: Steps 7 - 8 – Switch valve to LOAD position and fill loop with Inject Volume



**Figure 62.** Fast: Step 9 – Switch valve to INJECT position and start chromatographic process



**Figure 63.** Fast: Step 10 – Dispense rear sample volume and air segment to waste



## 5 Using Dynamic Load and Wash (DLW)

DLW Cycle Step-by-Step

**Figure 64.** Fast: Steps 11 - 12 – Clean valve with Wash Solvent 2




Figure 65. Fast: Steps 13 - 14 – Clean valve with Wash Solvent 1




Figure 66. Fast: Step 15 – Wash the syringe needle with Wash Solvent 1





Figure 67. Fast: End cycle





## Sample Trays

This appendix describes the nomenclature for the stack number, tray location number, and sample start location number.

**Note** You must use the handheld Control Terminal to specify the stack, tray location, and sample start numbers.

For more information about the trays, refer to the *Accela Open Autosampler Hardware Manual*.


### Contents

- [Tray Types](#)
- [Tray Position](#)
- [Sample Position](#)


## Tray Types

Figure 68 through Figure 71 show the tray types for use in the Accela Open Autosampler.


**Figure 68.** Tray type VT54



**Figure 69.** Tray type MT96




**Figure 70.** Tray type DW96



96-position, deep well micro plate (12 × 8)

**Figure 71.** Tray type MT384



384-position, high-density shallow micro plate  
(24 × 16)

## Tray Position


The nomenclature for the stack, tray, and sample positions is as follows:

CStk1–*TT*:*SS*

Where:

- CStk1 refers to the first stack cooler, which contains multiple tray drawers. (The designation for a second stack cooler is CStk2.)
- *TT* is the tray position within the overall stack. [Figure 72](#) shows the six tray positions, numbered 01 through 06, for a three-drawer stack.

**Figure 72.** Tray positions in a stack cooler with three drawers



- *SS* is the sample position within the designated tray. For examples, see [Table 17](#) and the [formula](#) that follows the table.

## Sample Position

For the 54-position tray, use the marked numbers on the tray to specify the sample start location. For the 96-position well plates, use the position numbers listed in [Table 17](#). For the 384-position well plate, use the [formula](#) that follows the table to calculate the sample start position number.

**Table 17.** Position numbers for the 96-position (12 × 8) micro plate

<b>A-B</b>	<b>No.</b>	<b>B-C</b>	<b>No.</b>	<b>C-D</b>	<b>No.</b>	<b>E-F</b>	<b>No.</b>	<b>F-G</b>	<b>No.</b>	<b>G-H</b>	<b>No.</b>
A1	1	B5	17	C9	33	E1	49	F5	65	G9	81
A2	2	B6	18	C10	34	E2	50	F6	66	G10	82
A3	3	B7	19	C11	35	E3	51	F7	67	G11	83
A4	4	B8	20	C12	36	E4	52	F8	68	G12	84
A5	5	B9	21	D1	37	E5	53	F9	69	H1	85
A6	6	B10	22	D2	38	E6	54	F10	70	H2	86
A7	7	B11	23	D3	39	E7	55	F11	71	H3	87
A8	8	B12	24	D4	40	E8	56	F12	72	H4	88
A9	9	C1	25	D5	41	E9	57	G1	73	H5	89
A10	10	C2	26	D6	42	E10	58	G2	74	H6	90
A11	11	C3	27	D7	43	E11	59	G3	75	H7	91
A12	12	C4	28	D8	44	E12	60	G4	76	H8	92
B1	13	C5	29	D9	45	F1	61	G5	77	H9	93
B2	14	C6	30	D10	46	F2	62	G6	78	H10	94
B3	15	C7	31	D11	47	F3	63	G7	79	H11	95
B4	16	C8	32	D12	48	F4	64	G8	80	H12	96

Use the following formula to calculate the sample position number:

$$N = ([Row - 1] \times W) + Col$$

Where:

- *N* is the calculated sample position number in the tray.
- *Row* is the sample's row converted to a number.  
A equals 1, B equals 2, C equals 3, and so on.
- *W* is the total number of columns on the tray:
  - Equals 12 for the 96-position tray
  - Equals 24 for the 384-position tray
- *Col* is the sample's column number:
  - Equals 1 to 12 for the 96-position tray
  - Equals 1 to 24 for the 384-position tray

The following example is for the sample position C22 in a 384-position micro plate:

$$N = ([3 - 1] \times 24) + 22 = 48 + 22 = 70$$



# Index

## A

- Accela Open Autosampler
  - adapter cable for detector 8
  - interconnect cable 6
- autosampler control terminal, using 19

## C

- cables
  - Accela Open Autosampler system interconnect cable 6
  - Accela Pump adapter cable 12
  - DB9 to DB9 serial cable 4
  - detector adapter cable 8
- compliance
  - FCC iii
  - regulatory iii
- compliance, WEEE v
- computer requirements xi
- contacting us xiii
- control terminal attached to the autosampler 19
- cycle composer macros 57
- cycle composer macros or ICC cycles, installing 57
- Cycle Editor 57

## D

- delay time 26
- device drivers xi
- Direct Control Configuration page 52
- Direct Control dialog box, opening 47
- DLW
  - characterizing features 55
  - cycle 68
  - cycle for fast injection 90
  - cycle for standard injection 68
  - fast injection cycle 63
  - macro definitions 57
  - operating 64
  - overview 55

- priming solvent lines 59
- standard injection cycle 60
- DLW Actuator/Solenoid 66
- DLW Pumps 66
- documentation
  - accessing x
  - additional x
  - online xiv
  - survey xiv

## E

- electromagnetic compatibility iii
- EMC compliance iii
- equation, sample position number in tray 108
- error message, server failure xi

## F

- FCC compliance iii
- file types
  - instrument methods (.meth) 26
  - macro (.pma) 41
  - templates (.pme) 32

## I

- ICC cycles 57
- Injected screen on hand-held controller 21
- Inlet Direct Control, opening 48
- Instrument Configuration window 15
- Instrument Setup window, opening 23

## J

- Job Queue screen 20

## L

level surface, importance of 2  
lithium battery 1

## M

Macro Editor  
  enabling 17  
  opening 39  
  page parameters 39  
Macro, Instrument Setup 28  
macros  
  creating 39  
  cycle composer 57  
  importing 17  
  list of standard 39  
mass spectrometry applications ix  
Menu screen 20  
method parameters, specifying 23  
method templates, list of 24  
MSQ Plus MS  
  Accela Open Autosampler system interconnect cable  
  connection 6, 9

## N

nomenclature for stacks, trays, and samples 103  
numeric variables 44

## O

object variables 44  
ordering information 2

## P

parameters  
  Direct Control  
  Method page  
  Description 49  
  Device Status 49  
  Injections Per Vial 50  
  Injectors 50  
  Method 49  
  Syringe 49  
  Trays 50  
  Vial Index 49  
  Volume 50  
Tray Configuration page  
  Tray 51  
  Tray Holder 51  
  Type 51

  Vials in X Direction 51  
  Vials in Y Direction 51  
Instrument Configuration 18  
  Auto Select 18  
  Browse 19  
  Check 18  
  Com Port 18  
  Enable 18  
  Import 19  
  Import Path 19  
  Overwrite Existing Macros 19  
  Stack 19  
  Tray Name 19  
Instrument Setup  
  Browse 27  
  Delay Time 28  
  Macro Sequence 27  
  Method Folder 27  
  New Template 27  
  Syringe 27  
  Template 27  
  Template Description 27  
  Variables 28  
  Volume 28

Macro Editor page  
  Add 41  
  Command 42  
  Command Description 43  
  Command Sequence 42  
  Delete 41–42  
  Edit 41  
  Insert Above 42  
  Insert Below 42  
  Macro 41  
  Macro Description 41  
  Macro Variables 41  
  Parameters 43  
  Replace 42

macros  
  Fast Injection Accela Open 63  
  Priming Accela Open AS\_Rev01 59  
  Stand Injection Accela Open 60  
Run Auto Sampler Method dialog box  
  Injectors 36  
  Trays 36  
  Vial Index 36  
  Volume 36  
Template Editor page 31  
Variable Definition dialog box  
  A Numeric Value 45  
  An Object 45  
  Default Value 45  
  Lower Limit 45  
  Upper Limit 45

- Variable Name 45
- Variable Units 45
- part numbers 2
- priming solvent lines 59
- pulse time of the out signals, setting 22
- pumps, self-priming 56

## R

- Recommended Injection Volume box 26
- regulatory compliance iii
- release notes xi
- requirements
  - computer xi
  - software xi
- Run Auto Sampler Method dialog box parameters 34

## S

- safety information 1
- safety standards iii
- sample trays 103
- self-priming micro pumps 56
- server failure error message xi
- software requirements xi
- solvent lines, priming 59, 65
- stack, cooler 106
- survey link xiv
- syringe size, specifying 26
- system requirements xi

## T

- Template Editor
  - enabling 17
  - opening 31
- Template Editor page parameters 31
- templates
  - editing 31
  - Single Step LC Injection 24
  - testing 34
- Thermo PAL driver for Xcalibur 57
- Tray Configuration page 50
- tray type, selecting with the PAL Control Terminal 20
- Trays 36
- trays
  - position in stack 106
  - sample position 107
  - type
    - DW96 105
    - MT384 105

- MT96 104
- VT54 104
- Trays screen 20
- Tune window, making single injections from 48

## V

- variables, defining 44

## W

- Waiting for contact closure status for LC pump 21
- WEEE compliance v

